



Package Content

1. CardiAid CT0207
2. CardiAid Protection Bag - CT0207P
3. Pre-connected Electrodes - CT0207E
4. Emergency Kit - CT0207EK

- Plastic Emergency Bag
- Mask
- Razor
- Gloves
- Scissors
- Alcohol Pad

CARDIAID CT0207 - TECHNICAL SPECIFICATIONS

Defibrillation System

- Operating Mode :** Automated (one button operation)
- Wave Form :** Biphasic, current-limited
- Shock Energy :** Low energy (max. 181J) or high energy (max.237J), programmable, energy adjusted to patient impedance.
- Shock Sequence :** Constant or escalating, programmable (factory setting)
- Charging Time :** Typically < 10 seconds*, maximum < 15 seconds*

ECG Analysis System

- Analysis Time :** Typically < 10 seconds*
- Derivation :** II
- Asystole Threshold :** < 0,2 mV
- Specificity NSR/Asystole :** > % 95**
- Sensitivity VF/VT :** > % 90**
- Reaction to Implanted**
 - Pacemaker :**Normal cardiac pacemaker rhythms are not detected as being shockable
- Impedance Measurement :** Checks electrodes contacts.
- Movement Detection :** Check the signal quality; acoustic warning at patient movement

Operation

- Operating Elements:**
 - i) Automatic switch-on when cover is opened
 - ii) Flashing shock key (one-button operation)
 - iii) Info key
- Info Mode :** Announcement of the elapsed time and number of shocks since device started, when info-button pushed
- Display Elements :** Self-explanatory illuminated symbols (traffic light principle)

Defibrillation Electrodes

- Range of Application :** CardiAid can be used for adult or young patients over 20 kgs.
- Delivery State :** Disposable self-adhesive electrodes ready for use, sealed and packed with connector outside of package
- Polarization :** Not polarized
- Electrode Surface Area :** Each 125 cm²
- Cable Length :** 130 cm
- Shelf Life :** 30 months from date of Manufacture

Data Management

- Event Documentation :** Automatic recording of ECG and event data in internal memory up to 4 sessions with a total maximum duration of 4 hours
- Data Transfer :** Bluetooth (only for service)
- Event Review and**
 - Device Programming :** Via bluetooth connection to PC with Software CardiAid View (only for service)

Self-Test

- Schedule :** Automatic daily, monthly and when device is switched-on
- Timing :** Programmable
- Scope :** Battery, internal electronics, software, charging

Energy Supply

- Type :** Alkaline
- Shock Capacity :** Typically 210 shocks*, min. 100 shocks*
- Monitoring Capacity :** Up to 20 hours
- Stand-by Period :** 2 years

Environmental Conditions

- Storage and Transport:**
 - i) 0°C to +50°C
 - ii) -20°C to +60°C (limited to max. 2 weeks, without battery and electrodes)
- Operation :** 0°C to +50°C
- Relative Humidity :** 0 - 95 %

Standards

- Norms :** EN60601-2-4, EN60601-1
- Resuscitation Protocol :** ERC, ILCOR 2005
- EMC :** EN60601-1-2, EN55011, EN61000-4-3, EN61000-4-8
- Protection Class :** IEC 529, IPX4 (spray water safety))

Dimensions and Weight

- Width :** 301 mm
- Height :** 304 mm (including handle)
- Depth :** 112 mm
- Weight :** 3,1 kgs (including battery and electrodes)
- Device Class :** II B (93/42/EEC)

* Data are valid at 20°C with new, fully-charged batteries. Values can vary and are dependent upon storage and environmental conditions, frequency of use, pre-configured settings, and the shelf life of the product.

** The algorithm and the board of CardiAid AED CT0207 offers the optimum combination of sensitivity and specificity.

1) Published in "Automated Analysis of electrical signals of the human body for detecting of lifethreatening cardiac abnormalities", Igor Tchoudovski, Mensch and Buch Verlag, ISBN 3-89820-984-9, Page 141

2) The New VF/VT algorithm and board is used in CardiAid AED CT0207.

CardiaTech, at its own discretion, reserves the right to make changes pertaining to the technical features and the accessories of the device.

CardiaTech, under no condition, shall be liable for any fault of spelling and/or act of writing in any and all instructive and informational data, documents, images given above and any direct, indirect, accidental, secondary or penal loss or damages may arise from them.